

Resistor-Temperature Conversion

Tesla cryogenic system local application

Wed, May 11, 1994

The TRTD local application supports needs of low temperature measurements for the Tesla superconducting RF laboratory in LAB-2 in the Fermilab village. Four signals are demultiplexed into 16, thermal emf's measured and used to get the resistance of the cold resistors, which are of two types, platinum and carbon. From the resistance, calibration tables are used with linear interpolation to get the temperature in degrees Kelvin. This note describes the methods used to do this.

Parameters layout

The available parameters for TRTD as seen via "Page E" are as follows:

```
E LOC APPL PARAMS 05/11/94 0838
NODE<056B>  NTRY< 5>  HELP<0576>
NAME=TRTD  CNTR=02AB
TITL"TEMPERATURE RESISTORS  "
SVAR=00046044
ENABLE  B<00C2> TRTD ENABLE
MPXDATA C<0120> T1MUX      V
V+,ETC  C<0200> HM1THV     V
R REF   C<01F0> RREFPT     OHM
MPXSEL  B<0128> RTD MPX SELECT 0
CURRDIR B<0120> RTD PT CURR DIR
EMF PER  <0258>
```

The usual enable Bit# argument is followed by MPXDATA, the base channel# of the four multiplexed signals. V+,ETC is the base channel# of the sixteen V+, V-, emf, resistance, and temperature values, altogether 5 sequential "arrays" of 16 channels. The R REF reference resistor channel# specifies the Pt reference resistor. The next channel# holds the Carbon reference resistor value. MPXSEL is the base Bit# of the digital byte used for selecting the four 1-of-4 multiplexer values. The program sequences this byte through values of 00,AA,55,FF on each four successive cycles. The CURRDIR parameter is the base Bit# of the two bits used for controlling the current direction through the Pt and C resistors, respectively. The EMF PER parameter is the number of cycles between measurements of the V- and emf values, done by reversing the current direction for four cycles to collect these data.

Demultiplexing

Signals from the RTD system are very small, so they are amplified by four-input multiplexed amplifiers. Two-bit select codes are wired to each amplifier. The program operates all four two-bit select codes at once, as they share the same byte of digital output. Data is read from one selection each cycle, which for the LAB-2 facility is 10 Hz. In this way, all 16 channels of data that result are updated every 400 ms.

Thermal emf's

To measure the thermal emf's that result from the conductor leads that connect from room temperature to cold temperatures, it is necessary to reverse the current. Every so often, perhaps every few minutes, the current is reversed and the data collected for the negative direction of current flow. After taking four cycles to accomplish this, the current is returned to the usual positive direction. The formula for the emf is as follows:

$$\text{emf} = (V^+ + V^-) / 2$$

These emf values are updated every time a new V- is measured with the current direction negative.

Resistances

The resistors used are of two types, platinum and carbon. The single Pt resistor calibration table is fairly linear down to liquid nitrogen temperatures, so linear interpolation can be used with a set of 10 calibration points only. These Pt calibration points are part of the program source and are as follows:

```
kelvinPt[1]:= 14.0;    ohmsPt[1]:= 1.797;
kelvinPt[2]:= 20.0;    ohmsPt[2]:= 2.147;
kelvinPt[3]:= 30.0;    ohmsPt[3]:= 3.508;
kelvinPt[4]:= 40.0;    ohmsPt[4]:= 5.938;
kelvinPt[5]:= 50.0;    ohmsPt[5]:= 9.228;
kelvinPt[6]:= 70.0;    ohmsPt[6]:= 17.128;
kelvinPt[7]:= 100.0;   ohmsPt[7]:= 29.987;
kelvinPt[8]:= 150.0;   ohmsPt[8]:= 50.815;
kelvinPt[9]:= 200.0;   ohmsPt[9]:= 71.073;
kelvinPt[10]:= 300.0;  ohmsPt[10]:= 110.452;
```

The formula for resistance is as follows:

$$R = ((V^+ - \text{emf}) / V^+_{\text{ref}}) * R_{\text{ref}}$$

Here V^+_{ref} is the voltage across the reference resistor (Pt or C) and R_{ref} the reference resistor value itself. The reference values are constants in the program. The values at the time of this writing are $R_{\text{ref}} = 99.82$ (Pt) and 1000.0 (C).

The carbon resistors calibration curve is very nonlinear, so more calibration points are used. Linear interpolation is done with 20 calibration points between the *log* of the temperature and the *log* of the resistance. In order to collect these calibration data for use by the TRTD program, the calibration data points were entered into Excel by editing the original full calibration data set text file. In the spreadsheet, it was easy to calculate an average of the 5 calibrated carbon resistors for use with those carbon resistors that have no measured calibration data. In addition, logs were calculated for the temperature and resistance calibration data points. These log values were output from Excel via a text file, which was then edited into an MPW assembly source file containing "DC.S" data directives. This source file was then assembled and the resulting data downloaded into the station that will run the TRTD local application. Upon initialization, the program requests this local data file and places the data into its calibration data arrays for use in converting the carbon resistances into temperatures. Part of the calibration data entered into Excel was this:

	T2	T3	T4	T5	T6	TAvg
940214						
293.477	107.54	106.33	107.70	108.29	101.33	106.24
77.828	127.02	125.42	126.61	127.68	123.28	126.00
56.171	139.35	137.44	138.79	139.94	135.31	138.17
46.159	148.57	146.50	147.92	149.07	144.34	147.28
...
4.686	823.95	794.60	815.01	812.69	789.51	807.15
4.254	960.07	923.74	948.93	945.27	918.79	939.36
3.849	1137.84	1092.06	1123.90	1118.58	1087.68	1112.01
2.597	2875.26	2729.90	2840.76	2798.62	2734.71	2795.85
2.083	5196.52	4879.04	5103.97	5036.64	4916.87	5026.61
1.813	7783.32	7266.99	7634.61	7521.43	7368.14	7514.90
1.527	13827.37	12689.70	13608.76	13288.20	13037.05	13290.22
1.257	28765.90	26017.65	28125.64	27273.40	26708.60	27378.24

The calculated logs of the calibration point data are as follows:

Log ₁₀ (temp)	Log ₁₀ (T2)	Log ₁₀ (T3)	Log ₁₀ (T4)	Log ₁₀ (T5)	Log ₁₀ (T6)	Log ₁₀ (TAvg)
2.4676	2.0316	2.0267	2.0322	2.0346	2.0057	2.0263
1.8911	2.1039	2.0984	2.1025	2.1061	2.0909	2.1004
1.7495	2.1441	2.1381	2.1424	2.1459	2.1313	2.1404
1.6643	2.1719	2.1658	2.1700	2.1734	2.1594	2.1681
...
0.6708	2.9159	2.9001	2.9112	2.9099	2.8974	2.9070
0.6288	2.9823	2.9655	2.9772	2.9756	2.9632	2.9728
0.5853	3.0561	3.0382	3.0507	3.0487	3.0365	3.0461
0.4145	3.4587	3.4361	3.4534	3.4469	3.4369	3.4465
0.3187	3.7157	3.6883	3.7079	3.7021	3.6917	3.7013
0.2584	3.8912	3.8614	3.8828	3.8763	3.8674	3.8759
0.1838	4.1407	4.1035	4.1338	4.1235	4.1152	4.1235
0.0993	4.4589	4.4153	4.4491	4.4357	4.4267	4.4374

These data were then edited as data declaration statements into an MPW assembly source file, assembled and downloaded as a data file of 32-bit floating point values called DATATRTD into the local station along with the LOOPTRTD local application. When the program is first initialized, it reads the data file for use during Carbon resistor temperature conversion.

Internal details

Four channels are read each cycle and converted as necessary. Of the 16 channels read, they are of the following "resistor types":

Signal#	Mpx	Resistor type
1-4.	00	Pt, Pt, Pt, Pt
5-8.	01	CAvg, CAvg, CAvg, CAvg
9-12.	10	Pt, Pt, none, none
13-16.	11	C6, C2, C3, C4

The resistor type specifies which reference resistor is used as well as which calibration data to use to perform the temperature conversion.

During program operation, a data structure allocated at initialization keeps track of the context. Using Page E, one can find the location of this dynamically assigned data structure. Its structure at the time of this writing, expressed in the format produced by the memory dump page application, at a time when the signal values may not have been valid, is as follows:

```

M MEMORY DUMP      05/12/94 1056
056B:03D716 0420 0000 0000 0007   BlkSize, , ,BlkType
056B:03D71E 056B 0004 0000 0000   Node#, state#, wait, debug#
056B:03D726 0000 0085 0200 0210   currDir, dirWait, V+Chan, V-Chan
056B:03D72E 0220 0230 0240 1004   emfChan, resistChan, tempChan, dbI
056B:03D736 0000 00AA 0055 00FF   mpx, mpx, mpx, mpx
056B:03D73E 0000 0CC7 0000 7FFF   resRef(Pt), resRef(C)
056B:03D746 467F 8C00 44CC CB33   refConst(Pt), refConst(C)
056B:03D74E 7368 7480 03C3 0013   mpxChan(1-4)

056B:03D756 0101 0101 0707 0707   resTypes(1-16)
056B:03D75E 0101 0000 0602 0301

```

056B:03D766	0103	0102	0202	0102	deltaTimes(1-16)
056B:03D76E	0102	0102	0102	0202	
056B:03D776	73EF	73CD	73AC	73C6	V+(1-16)
056B:03D77E	7374	7493	7398	73FD	
056B:03D786	7393	7400	03C9	001C	
056B:03D78E	73AC	73FE	73CA	73D4	
056B:03D796	7387	73BB	73AC	73BD	V-(1-16)
056B:03D79E	738D	73E7	738B	73B9	
056B:03D7A6	7374	7402	03C2	0021	
056B:03D7AE	73A7	73F2	7249	738D	
056B:03D7B6	7382	73D1	73B6	73D3	emf(1-16)
056B:03D7BE	7391	73F2	7386	73F9	
056B:03D7C6	7380	73EE	03B5	001B	
056B:03D7CE	73BC	73E8	72EB	7434	
056B:03D7D6	016B	0604	FDB4	FEC8	resist(1-16)
056B:03D7DE	08BA	E333	099A	EF6B	
056B:03D7E6	0140	023D	0000	0000	
056B:03D7EE	F71D	0AAB	3FA4	FD91	
056B:03D7F6	0540	0CE9	0000	0000	kelvin(1-16)
056B:03D7FE	00E6	0042	00DC	0042	
056B:03D806	0445	07A9	0000	0000	
056B:03D80E	0042	00D5	006C	0000	
056B:03D816	3FE6	0419	4009	6873	Pt calibration points
056B:03D81E	4060	8312	40BE	0419	
056B:03D826	4113	A5E3	4189	0625	
056B:03D82E	41EF	E560	424B	428F	
056B:03D836	428E	2560	42DC	E76D	
056B:03D83E	4160	0000	41A0	0000	
056B:03D846	41F0	0000	4220	0000	
056B:03D84E	4248	0000	428C	0000	
056B:03D856	42C8	0000	4316	0000	
056B:03D85E	4348	0000	4396	0000	
056B:03D866	0000	0000	0000	0000	C calibration points
056B:03D86E	0000	0000	0000	0000	
...	
056B:03D906	4002	05BC	4006	A64C	Log10(Temp)
056B:03D90E	4009	38EF	400B	0069	
056B:03D916	400D	7F63	4011	AEE6	
056B:03D91E	4014	F41F	4019	B08A	
056B:03D926	401E	D5D0	4026	8241	
056B:03D92E	402E	0C4A	4037	DD98	
056B:03D936	403A	9E1B	403E	DE01	
056B:03D93E	4043	9724	405D	5B57	
056B:03D946	406D	CE07	4079	096C	
056B:03D94E	4084	809D	408E	AF4F	
056B:03D956	4001	B574	4006	4C30	Log10(T2)
056B:03D95E	4008	D6A1	400A	9C78	
056B:03D966	400D	119D	4011	35A8	
056B:03D96E	4014	6DC6	4019	1EB8	
056B:03D976	401E	339C	4025	C5D6	

056B:03D97E	402D	374C	4036	E48F	
056B:03D986	4039	9B3D	403D	CAC1	
056B:03D98E	4042	71DE	405B	E910	
056B:03D996	406C	0D1B	4077	212D	
056B:03D99E	4083	4FDF	408D	4A23	
056B:03D9A6	4002	0F91	4006	8F5C	Log10(T3)
056B:03D9AE	4009	1D15	400A	E148	
056B:03D9B6	400D	6388	4011	8ADB	
056B:03D9BE	4014	CB29	4019	85F0	
056B:03D9C6	401E	A993	4026	4A8C	
056B:03D9CE	402D	CAC1	4037	923A	
056B:03D9D6	403A	511A	403E	8A72	
056B:03D9DE	4043	3EAB	405D	0481	
056B:03D9E6	406D	4E3C	4078	7FCC	
056B:03D9EE	4084	4817	408E	5F07	
056B:03D9F6	4002	36E3	4006	CA58	Log10(T4)
056B:03D9FE	4009	566D	400B	18FC	
056B:03DA06	400D	94AF	4011	B08A	
056B:03DA0E	4014	F27C	4019	A1CB	
056B:03DA16	401E	B9F5	4026	511A	
056B:03DA1E	402D	C91D	4037	81D8	
056B:03DA26	403A	3BCD	403E	703B	
056B:03DA2E	4043	1DE7	405C	9A02	
056B:03DA36	406C	EF35	4078	154D	
056B:03DA3E	4083	F3B6	408D	F141	
056B:03DA46	4000	5D64	4005	D14E	Log10(T5)
056B:03DA4E	4008	6738	400A	339C	
056B:03DA56	400C	AF4F	4010	D35B	
056B:03DA5E	4014	13A9	4018	C2F8	
056B:03DA66	401D	DB23	4025	7C1C	
056B:03DA6E	402C	F5C3	4036	B368	
056B:03DA76	4039	6F00	403D	A512	
056B:03DA7E	4042	5604	405B	F62B	
056B:03DA86	406C	44D0	4077	837B	
056B:03DA8E	4083	AFB8	408D	A787	
056B:03DA96	4001	AEE6	4006	6CF4	Log10(T6)
056B:03DA9E	4008	FC50	400A	C227	
056B:03DAA6	400D	3F7D	4011	652C	
056B:03DAAE	4014	A3D7	4019	59B4	
056B:03DAB6	401E	7525	4026	147B	
056B:03DABE	402D	8FC5	4037	4F0E	
056B:03DAC6	403A	0C4A	403E	425B	
056B:03DACE	4042	F34D	405C	9375	
056B:03DAD6	406C	E219	4078	0EBF	
056B:03DADE	4083	F3B6	408D	FF2E	
056B:03DAE6	401D	ED29	3FF2	0F91	Log10(TAvg)
056B:03DAEE	3FDF	EF9E	3FD5	07C8	
056B:03DAF6	3FC7	65FE	3FB3	CD36	
056B:03DAFE	3FA6	D289	3F96	CBFB	
056B:03DB06	3F88	4817	3F6D	5CFB	

```
056B:03DB0E 3F51 4120 3F33 126F
```

```
056B:03DB16 3F2B B98C 3F20 F909
```

```
056B:03DB1E 3F15 D639 3ED4 3958
```

```
056B:03DB26 3EA3 2CA5 3E84 4D01
```


```
056B:03DB2E 3E3C 3611 3DCB 5DCC
```

The diagnostic `deltaTimes` byte array shows elapsed times per 10 Hz cycle in units of 0.5 ms. (A value of 02 means 1 ms.) Since the program operates on 4 signals per cycle, these times exhibit a periodicity of four cycles.

All floating point values are in IEEE standard 32-bit single precision.

Debug mode

A special debug mode is available for checking the ohms to temperature conversion for the Carbon resistors. To use it, find the above structure using Page E. Then set the `debug#`, which is initialized to zero, to the index value 1–16 of the channel whose conversion is to be checked. When this word is nonzero in this range, the program stops doing its usual work and concentrates on only the given channel index. It accepts the resistance value from the resistance channel, rather than computing a new one based upon the current demultiplexed data readings. This allows a value to be set into a resistance channel and the temperature channel observed resulting from this conversion. One can make a calibration plot of temperature versus resistance in this way using, for example, the Parameter Page on the Macintosh written by Bob Peters. Here is an example of making such a calibration curve for the T6 signal:



Logarithms, etc

Because of the extremely nonlinear characteristic of the resistance of the Carbon resistors relative to temperature, it was necessary to work with logarithms of these parameters. With linear interpolation on the *logarithms*, reasonably accurate results can be obtained. But the IRMs use a 68040 cpu, which does not have support for logarithms and powers on-chip, although it does have support for the basic add, subtract, multiply and divide. So a routine was needed to compute logarithms and powers. In a book called "Approximations for Digital Computers" by Cecil Hastings, Jr., Princeton University Press, 1955, can be found suitable formulae for these and other functions to several degrees of accuracy. The formula for $\text{Log}_{10}(x)$ that was used in TRTD has an error function < 0.000002 over the range $1 \leq x \leq 10$. Its form is as follows:

$$\text{Log}_{10}(x) = 0.5 + c_1*f + c_3*f^3 + c_5*f^5 + c_7*f^7,$$

where $f = (x - k)/(x + k)$, and $k = \sqrt{10} = 3.162278$. The set of coefficients is:

$$c_1 = 0.86855434; \quad c_3 = 0.29115068; \quad c_5 = 0.15361371; \quad c_7 = 0.21139497.$$

The formula for 10^x that was used has an error function < 0.000003 over the range $0 \leq x \leq 1$. Its form is as follows:

$$10^x = (1 + a_1x + a_2x^2 + a_3x^3 + a_4x^4 + a_5x^5)^2$$

The set of coefficients is:

$$a_1 = 1.15138424; \quad a_2 = 0.66130851; \quad a_3 = 0.26130650; \quad a_4 = 0.05890681; \quad a_5 = 0.02936622;$$

The TRTD program has about 600 lines of Pascal and compiles into about 4K bytes.